

Reports for Monitoring AUC Utilization

The following reports can be used to monitor the utilization of AUC checking for advanced imaging procedures to ensure compliance with the PAMA Imaging Appropriate Use Criteria requirement.

AUC Utilization Reports

The AUC utilization reports can be run by category and date range to generate a list of orders where AUC values were captured. Links to the Report Writer archive are provided in the table below. From the Report Writer archive links provided you can download the report and then upload into your MEDITECH system.

Technology	Report Link
Expense, 6.x, & C/S	ID#02305
\$T RAD	ID#02276
MG ITS	ID#02277

Note: After you download this file, you upload it into your MEDITECH system. Do NOT open or unzip this file using any sort of unzip utility, refer to KB article [#57720](#) for detailed instructions.

Deficiency Reports

AUC Deficiency Reports can be run by Facility, Registration Type, Financial Class, and Category for AUC-enabled procedures to identify orders where AUC values were not captured.

Order Management > Clinical > Reports > Imaging AUC Deficiency Report

AUC Deficiency Reports are available for every platform and can be accessed using the paths below.

Platform	Report Path
MG	Ancillary Dep > Imaging AUC Deficiencies
C/S	Care Area > Order Reports > Imaging AUC Deficiency Report
6x/Expense	Clinical > Reports > Ordering > Imaging AUC Deficiency Report