


# MEDITECH's

## Anesthesia Interface Suite



### Solutions for Safe Anesthesia Management

MEDITECH has joined forces with Spectrum Medical's Anesthesia Information Management System (AIMS) to create a suite of anesthesia interfaces which provide care teams with a complete clinical and medication record to support patient safety and improved outcomes. Additionally, workflow automation, increased charge capture, and time-saving processes can provide extensive financial savings.



(The image above depicts a sample of data exchange between MEDITECH Expance and Spectrum Medical's AIMS)

### What's in the Anesthesia Interface Suite?

- Outbound from MEDITECH: Medications (Administered and Home), ADT, Allergies, Labs, Scheduling, Vital Signs.
  - Inbound to MEDITECH: Medications (Administered\* and Home), Allergies, Anesthesia Reports, Billing, Intake & Output Data, Problems, Vital Signs.
- \*Administered medications inbound is available with Expance only.*
- Direct launching point access between systems.

### About Spectrum Medical



**Spectrum Medical** supplies anesthesia and perfusion solutions that are specifically designed and configured for the clinical specialty. Solutions include electronic medical record systems, comprehensive medical device connectivity, hardware, telemedicine, and reporting capabilities that are designed to complement MEDITECH. The new suite of anesthesia interfaces provides two-way communication of clinical and patient information to help ensure patient safety, quality of care, and efficiency of care.

### Learn More

For more information regarding deployment of these anesthesia interfaces, pricing, and product details, please contact your MEDITECH marketing consultant at: (781) 821-3000.